

AVS-N SERIES - INSTALLATION GUIDE

Information to consider before installing your RAYPA autoclave.

INDEX

AVS-8 MODEL

Electrical connection	Page 2
Connections graph	Page 2
Components included	Page 3
Drainage connection	Page 4
Water supply	Page 4
Dimensions to consider	Page 5
Environmental conditions	Page 5
Recommended maintenance	Page 6

AVS-12 MODEL

Electrical connection	Page 7
Connections graph	Page 7
Components included	Page 8
Drainage connection	Page 9
Water supply	Page 9
Dimensions to consider	Page 10
Environmental conditions	Page 11
Recommended maintenance	Page 11




TOP-LOADING BENCHTOP
AUTOCLAVES WITHOUT
DRYING
AVS-N SERIES

AES-8 MODEL


ELECTRICAL CONNECTION

The following table shows the plug configuration according to IEC and SCHUKO international standards. For customers who require other plugs and other electrical configurations, please contact our technical team.

MODEL	FREQUENCY	POWER	AMPERES / PHASE	VOLTAGE	CONNECTION
AES-8	50/60 Hz	1000 W	4 A	230 (1P+N+PE) V	16 A ①
AES-8-115V	50/60 Hz	1000 W	8 A	120 (1P+N+PE) V	16 A ①



CONNECTIONS GRAPH




TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-8 MODEL

COMPONENTS INCLUDED

In addition to the accessories chosen at the time of purchase of the autoclave the following components are included:



1 stainless steel protective cover with legs for the heating elements of Ø208 mm to place on the inner base of the sterilization chamber.



1 stainless steel wire basket of Ø220 x 150 mm.



1 heat resistant silicone hose of Ø8 x Ø14 mm and 1 m long for draining the sterilization chamber.

For:

3. STERILIZATION CHAMBER DRAIN OUTLET


TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-8 MODEL

DRAINAGE CONNECTION

To drain the water from the sterilization chamber, connect the supplied silicone hose to a drain or place a container (such as a tray) in the STERILIZATION CHAMBER DRAIN OUTLET **(3)**.

See "DRAIN HEIGHT of the sterilization chamber" on page 5.




The valve **(1)** must be rotated to the "DRAIN" position to drain the sterilization chamber.

WATER SUPPLY FOR STERILIZATION

Purified water is necessary for the operation of the autoclave and is added directly and manually into the sterilization chamber until it covers completely the heating elements cover (around 1.2 L).

The valve **(1)** must be in the "CLOSED" position when filling and during regular operation of the autoclave.



ECOPUR-500 accessory

The ECOPUR-500 accessory is an eco-efficient purifier designed to filter water directly from a water inlet, allowing for easy and economical filling of containers. This system is ideal as an alternative in situations where purified water networks are unavailable or when purchasing purified water containers is not preferred. With high-efficiency reverse osmosis technology, our ECOPUR-500 water purifier removes up to 99% of chlorine and chemical contaminants, ensuring high-quality water. Additionally, it features a quickconnect system for easy and secure installation, making this accessory a practical and sustainable solution for providing purified water anywhere.



IMPORTANT NOTE:

PURIFIED WATER

The water quality used to feed the autoclave must be free of contaminants and meet the following hardness and conductivity requirements:

- Hardness: ≤ 0.02 mmol/l
- Conductivity: between 5 μ S/cm and 15 μ S/cm

Multiple systems can be used to obtain water with these requirements: osmotized, demineralized, decalcified, distilled water, etc*.

*Note: Take into account that distilled water that is too pure (conductivity less than 5 μ S/cm) is not recommended as it may cause long term corrosion on stainless steel.

**TOP-LOADING BENCHTOP
AUTOCLAVES WITHOUT
DRYING**
AVS-N SERIES


AES-8 MODEL




DIMENSIONS TO CONSIDER FOR THE INSTALLATION OF YOUR AUTOCLAVE

For safety reasons, the distance between both sides of the autoclave and the wall or any other object must be 100 mm, and between the autoclave and the rear wall must be at least 200 mm.


MODEL	L LENGTH	D DEPTH	D1 DEPTH with maximum door opening	H HEIGHT	H1 HEIGHT with maximum door opening	LA x DA SUPPORT AREA	HD DRAINAGE HEIGHT of the sterilization chamber
AES-8	410 mm	355 mm	415 mm	430 mm	625 mm	272 x 272 mm	50 mm



WARNING:
Respect the recommended distances.



FRONTAL AND LATERAL VIEWS



**ENVIRONMENTAL
CONDITIONS**

This autoclave can operate under the following maximum conditions:

- Ambient temperature: 30 °C
- Humidity: 75%
- Altitude: 3,000 meters above sea level. Take into account that from 1,000 meters above sea level an adjustment of the purge parameter must be made; consult with our technical team.

TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-8 MODEL

RECOMMENDED MAINTENANCE

Autoclaves are like cars: they need regular maintenance for them to work properly, to ensure good condition, prevent deterioration of its components, and maximize their useful life. Frequent cleaning and regular maintenance are essential, as autoclaves work at high pressures and temperatures and are therefore subject to a high level of stress.

For clients that perform multiple cycles each day, we recommend carrying out maintenance and cleaning tasks more frequently.

DAILY MAINTENANCE

Clean the gasket using a clean cotton cloth with a mild vinegar solution (or a comparable product) to reduce the formation of lime deposits.

Wipe the external surfaces with a clean cotton cloth moistened with water and a neutral detergent. Ensure all surfaces are thoroughly dried afterwards.



WEEKLY MAINTENANCE

Clean the sterilization chamber and all the components with a clean cotton cloth dampened with water and neutral detergent. Dry them after.

Perform a complete drainage of the sterilization chamber and proceed to refill it with purified water.



ANNUAL MAINTENANCE

Conducting a technical inspection is essential to ensure consistent process safety over time. It is necessary to periodically verify the thermodynamic process parameters (pressure and temperature), ensuring they remain within the accepted limits.



TOP-LOADING BENCHTOP
AUTOCLAVES WITHOUT
DRYING
AVS-N SERIES

AES-12 MODEL


ELECTRICAL CONNECTION

The following table shows the plug configuration according to IEC and SCHUKO international standards. For customers who require other plugs and other electrical configurations, please contact our technical team.

MODEL	FREQUENCY	POWER	AMPERES / PHASE	TENSION	CONNECTION
AES-12	50/60 Hz	1000 W	4 A	230 (1P+N+PE) V	16 A ①
AES-12-115V	50/60 Hz	1000 W	8 A	120 (1P+N+PE) V	16 A ①



CONNECTIONS GRAPH




TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-12 MODEL

COMPONENTS INCLUDED

In addition to the accessories chosen at the time of purchase of the autoclave the following components are included:



1 stainless steel protective cover with legs for the heating elements of Ø245 mm to place on the inner base of the sterilization chamber.



1 stainless steel wire basket of Ø220 x 200 mm.



1 heat resistant silicone hose of Ø8 x Ø14 mm and 1 m long for draining the sterilization chamber.

For:

5. STERILIZATION CHAMBER DRAIN OUTLET


TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-12 MODEL

DRAINAGE CONNECTION

To drain the water from the sterilization chamber, connect the supplied silicone hose to a drain or place a container (such as a tray) in the STERILIZATION CHAMBER DRAIN OUTLET (5).

See "DRAIN HEIGHT of the sterilization chamber" on page 10.




The valve (2) must be rotated to the "DRAIN" position to drain the sterilization chamber.

WATER SUPPLY FOR STERILIZATION

Purified water is necessary for the operation of the autoclave and is added directly and manually into the sterilization chamber until it covers completely the heating elements cover (around 2 L).

The valve (2) must be in the "CLOSED" position when filling and during regular operation of the autoclave.



ECOPUR-500 accessory

The ECOPUR-500 accessory is an eco-efficient purifier designed to filter water directly from a water inlet, allowing for easy and economical filling of containers. This system is ideal as an alternative in situations where purified water networks are unavailable or when purchasing purified water containers is not preferred. With high-efficiency reverse osmosis technology, our ECOPUR-500 water purifier removes up to 99% of chlorine and chemical contaminants, ensuring high-quality water. Additionally, it features a quickconnect system for easy and secure installation, making this accessory a practical and sustainable solution for providing purified water anywhere.



IMPORTANT NOTE:

PURIFIED WATER

The water used to feed the autoclave must be free of contaminants and meet the following hardness and conductivity requirements:

- Hardness: ≤ 0.02 mmol/l
- Conductivity: between 5 μ S/cm and 15 μ S/cm

Multiple systems can be used to obtain water with these requirements: osmotized, demineralized, decalcified, distilled water, etc*.

*Note: Take into account that distilled water that is too pure (conductivity less than 5 μ S/cm) is not recommended as it may cause long term corrosion on stainless steel.

TOP-LOADING BENCHTOP
AUTOCLAVES WITHOUT
DRYING
AVS-N SERIES


AES-12 MODEL




DIMENSIONS TO CONSIDER FOR THE INSTALLATION OF YOUR AUTOCLAVE

For safety reasons, the distance between both sides of the autoclave and the wall or any other object must be 100 mm, and between the autoclave and the rear wall must be at least 200 mm.


MODEL	L LENGTH with closed door	L1 LENGTH with maximum door opening	D DEPTH	D1 DEPTH with maximum door opening	H HEIGHT	HL HEIGHT OF LOAD	LA x DA SUPPORT AREA	HD DRAINAGE HEIGHT of the sterilization chamber
AES-12	490 mm	780 mm	475 mm	600 mm	630 mm	400 mm	437 x 412 mm	50 mm




WARNING:
Respect the recommended
distances.



FRONTAL AND LATERAL VIEWS



PLAN VIEW



TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

AES-12 MODEL

ENVIRONMENTAL CONDITIONS

This equipment can operate under the following maximum conditions:

- Ambient temperature: 30°C
- Humidity: 75%
- Altitude: 3000 meters above sea level. Please note that from 1000 meters above sea level, an adjustment of the purge parameter must be made. Consult with RAYPA technical team.

RECOMMENDED MAINTENANCE

Autoclaves are like cars: they need regular maintenance for them to work properly, to ensure good condition, prevent deterioration of its components, and maximize their useful life. Frequent cleaning and regular maintenance are essential, as autoclaves work at high pressures and temperatures and are therefore subject to a high level of stress.

For clients that perform multiple cycles each day, we recommend carrying out maintenance and cleaning tasks more frequently.

DAILY MAINTENANCE

Clean the gasket using a clean cotton cloth with a mild vinegar solution (or a comparable product) to reduce the formation of lime deposits.

Wipe the external surfaces with a clean cotton cloth moistened with water and a neutral detergent. Ensure all surfaces are thoroughly dried afterwards.



WEEKLY MAINTENANCE

Clean the sterilization chamber and all the components with a clean cotton cloth dampened with water and neutral detergent. Dry them after.

Perform a complete drainage of the sterilization chamber and proceed to refill it with purified water.



ANNUAL MAINTENANCE

Conducting a technical inspection is essential to ensure consistent process safety over time. It is necessary to periodically verify the thermodynamic process parameters (pressure and temperature), ensuring they remain within the accepted limits.



TOP-LOADING BENCHTOP AUTOCLAVES WITHOUT DRYING AVS-N SERIES

TECHNICAL SUPPORT, ORIGINAL SPARE PARTS AND EXPERIENCED CONSULTING SERVICES

For an optimal operation of the autoclave, always use original spare parts and schedule periodically, with a specialized technician, the relevant maintenance tasks, such as calibration of the temperature and pressure sensors or replacement of the sterilization chamber door gasket.

Additionally, we provide a comprehensive range of services to ensure a satisfactory user experience throughout the entire lifespan of our products. These services include support and training programs, guided start-up and qualification services, preventive and corrective maintenance, periodic calibration, technical support and repairs, as well as consulting on special applications.

Should you encounter any issues, have questions, or require further information regarding maintenance of the AVS-N Series autoclave, please don't hesitate to reach out to our technical support service using the following contact details.



Technical support

<https://www.raypa.com/en/contact/>
+34 937 830 720 (Ext. 2109)

more info

 YouTube



CLICK!
ACCESS ALL
VIDEOS FROM
RAYPA

Discover more information about our products on our **Youtube Channel**



REV 09/2024